Tracheostom Product Catalogue 2025

Bringing the right products and care your way **Care**

European Contact Centre & Customer Service Awards

2023

SILVER WINNER

Contents

Accessories18Neckbands19Dressings21Swabs24Wipes and Creams26Shower Protectors and Tape28Lubricating jelly and irrigation fluid29Dressing packs30Gloves30	HMEs and Speaking Valves	10
Dressings21Swabs24Wipes and Creams26Shower Protectors and Tape28Lubricating jelly and irrigation fluid29Dressing packs30	Accessories	18
Protectors	Dressings Swabs Wipes and Creams Shower Protectors and Tape Lubricating jelly and irrigation fluid Dressing packs Gloves	21 24 26 28 29 30 30

Kids	36
HMEFs Neckbands	37 39
Jaw Rehabilitation	44





Welcome to Atos Care

At Atos Care, our mission is to help make life easier for people living with a neck stoma. We understand that having a laryngectomy or tracheostomy can be life changing, and we are here to support you every step of the way.

In this catalogue you will find information on the services available to you and many of the products we supply. If you can't find what you're looking for, please contact our helpful Customer Care team, who can help you find the products you need if they are available on prescription. We supply products made in-house by Atos (the Provox Life, Provox and Tracoe ranges) and products from a variety of other manufacturers, to give you the best choice and selection.

We are here to give you a helping hand, whatever you may need. If you have any questions or feedback, please do not hesitate to get in touch.

How to order

At Atos Care, we are committed to making life easier for people living with a neck stoma. We offer a selection of options for how you can order your supplies; ranging from a monthly service call to a fully automated monthly order.

A friendly monthly service call, email or text

Our helpful Customer Care team can contact you each month to help you place your order and answer any questions you may have. Contact us when you wish to order

If you prefer to contact us when you need to order, you can contact us using our convenient freephone phone line, or by text, email or through our website, if you prefer. Automated convenience order

We can set up a convenient repeat order for you, so that your supplies are automatically ordered and sent to you each month.

Discover our new wesbite www.atos-care.co.uk for an easy and convenient way to order products

Order process explained

- 1
- Place an order with our helpful Customer Care team.
- 2

3

- If your order contains new items, we may check suitability with your healthcare professional.
- We contact your GP to request the prescription for your items.*
- 4
- On receiving the prescription**, we send out your supplies directly to your home.

Please allow up to 14 working days from placing your order to receiving your supplies.

The Catalogue doesn't contain everything we can provide, so if there is anything you need but cannot see, then please get in touch with your Customer Care Representative.

^{*}This applies if your GP agrees to receive prescriptions from a third party (i.e. us). If this is not the case, you may be required to place the request directly with your GP surgery.

^{**}Our Electronic Prescriptions Service (EPS) means your GP can send the prescription directly through us through a secure system, so that you don't have to collect it and send it on to us.

Events

We run a range of events, both online and in-person. Our MyLife events are one of the ways we can help support you to breathe better, whatever you do.

We hold a range of helpful webinars, practical workshops and interactive experiences. Our MyLife events include helpful presentations about relevant topics, as well as a chance to ask questions and hear from other people who know and understand about tracheostomy.

Access to experts

Supported events can also include healthcare and other professionals who can share valuable information, tips and best practice that might be hard to find on your own. Past events have included presentations and speakers on themes such as:

- Looking after your lungs
- Travel
- Breathing exercises
- Dealing with the unexpected









Breathing, speaking and living well

Tracoe[®] is a tracheostomy brand that includes a range of high quality devices designed to meet high standards for comfort, reliability and durability. HMEs (Heat and Moisture Exchangers) help to restore the lost function of the nose, enabling those with a tracheostomy to breathe better and to speak again.



TrachPhone® HME

TrachPhone is a versatile and lightweight HME that provides good humidification for different situations.

- A multifunctional HME offering effective lung rehabilitation
- Speaking valve function, press the top lid to occlude the stoma for voicing
- Integrated oxygen supply port (4 mm)
- Integrated suctioning port that allows for tracheal suction without removing the HME

ltem No.	Supplier	Description	Pack size
7707	Atos Medical AB	TrachPhone HME	30



XtraCare[™] HME-F

XtraCare HME provides good humidification and effective filtration, for extra protection.

- Combines an HME with an electrostatic filter, for extra protection
- Protects against airborne particles, including viruses and bacteria, with a filtration efficiency of more than 99%*
- Can be used with Freevent O² Adaptor for administration of supplemental oxygen

Item No.	Supplier	Description	Pack size
7767	Atos Medical AB	XtraCare HME White	30



XtraCare[™] HME-F Blue

Blue edition is developed to make it easy to see when you are using XtraCare HME-F.

ltem No.	Supplier	Description	Pack size
7768	Atos Medical AB	XtraCare HME Blue	30



Freevent® HME

Heat and Moisture Exchanger to be used with DualCare Speaking Valve or DigiTop.

ltem No.	Supplier	Description	Pack size
7742	Atos Medical AB	Freevent HME 15 Regular	30



O² Adaptor

- Adaptor for administration of supplemental oxygen
- To be used with XtraCare

ltem No.	Supplier	Description	Pack size
7769	Atos Medical AB	O2 Adaptor	10



DualCare[™] Speaking Valve

Combining a speaking valve with the benefits of a Heat and Moisture Exchanger (HME), with a simple twist.

- Unique combination of a speaking valve with a Heat and Moisture Exchanger (HME)
- Change between speaking mode and HME mode with a simple twist
- Hands-free speech in the speaking mode
- Improves or maintains lung health in HME mode

Item No.	Supplier	Description	Pack size
7744	Atos Medical AB	DualCare Speaking Valve	1



DualCare[™] Speaking Valve Blue

Blue edition is developed to make it easy to see when you are using a DualCare Speaking Valve.

ltem No.	Supplier	Description	Pack size
7755	Atos Medical AB	DualCare Speaking Valve Blue	1



Trachie-Voice Speaking Valve

ltem No.	Supplier	Description	Pack size
TSV/100	Insight Medical	Trachie-Voice Speaking Valve	2



Passy Muir® Valve

Item No.	Supplier	Description	Pack size
TRPMV1002	Passy Muir	Passy Muir Valve 2000 (clear)	1
TRPMV1003	Passy Muir	Passy Muir Valve 2001 (purple)	1
TRPMV0002	Passy Muir	Passy Muir Valve 2001 (Aqua)	1





Accessories

Accessories - Neckbands

Freevent® Neckband One-piece with Twill Ties

- One-piece neckband with twill ties, designed for extra security
- Simply cut and adjust to the exact size needed
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort

ltem No.	Supplier	Description	Size	Pack size
1661	Atos Medical AB	Freevent Neckband - One-piece	Upto 550mm	10

Freevent® Neckband Two-piece with VELCRO® fastening

- High quality adjustable neckband with VELCRO[®] fastening designed to provide secure tube positioning
- Simple application with VELCRO® fastening, adjustable to different neck sizes
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort

ltem No.	Supplier	Description	Size	Pack size
1662	Atos Medical AB	Freevent Neckband - Two-piece	Upto 550mm	10







Trachi-Hold® Tube Holder

Item No.	Supplier	Description	Size	Pack size
TR ACC 0014	Sentient	Trachi-Hold Small Tube Holder	Upto 350mm	10
TR ACC 0012	Sentient	Trachi-Hold Large Tube Holder	Upto 470mm	10



Insight Tube Holder Two-piece

ltem No.	Supplier	Description	Size	Pack size
TH100	Insight Medical	Insight Tube Holder - Two-piece	Upto 470mm	10

Marpac[™] Tube Holder with VELCRO[®]/twill hybrid fastener



Item No.	Supplier	Description	Size	Pack size
201	CCMed	Marpac Tube Holder - Two-piece	Upto 406mm	10

Marpac[™] Tube Holder with Twill Ties



ltem No.	Supplier	Description	Size	Pack size
108D	CCMed	Marpac Tube Holder - Two-piece	Upto 530mm	25

Accessories - Dressings



Tracoe Purofoam

Tracoe Purofoam dressing is made of a highly absorbent open pore foam. The outer layer is covered with a beige coating which prevents leakage of secretions to the outside and reduces bacterial contamination.

- Allows for quick absorption of secretions
- · Non-stick material for easy removal and replacement
- Circular opening and zigzag slit enables easy replacement whilst wearing the tracheostomy tube

Item No.	Supplier	Description	Size	Pack size
Ref 958	Atos Medical AB	Tracoe Purofoam Large	90mm x 100mm	10
Ref 959	Atos Medical AB	Tracoe Purofoam Small	65mm x 65mm	10

Trachi-Dress

- Slimline cotton-type dressing, suitable for those with minimal secretions
- Fenestrated with non-adhesive backing
- Non-sterile dressing

ltem No.	Supplier	Description	Size	Pack size
TRDRE1001	Atos Medical AB	Trachi-Dress Dressing Small	60mm x 80mm	20
TRDRE1002	Atos Medical AB	Trachi-Dress Dressing Large	80mm x 100mm	20



Metalline Dressing

- Slimline cotton-type dressing, suitable for those with minimal secretions
- Aluminised backing designed to promote wound healing

ltem No.	Supplier	Description	Size	Pack size
23 094	Lohmann & Rauscher	Metalline Dressing	80mm x 90mm	50



ActivHeal

- Highly absorbent foam dressing, suitable for those with lots of secretions
- Fenestrated with non-adhesive backing

Item No.	Supplier	Description	Size	Pack size
10009118	Advanced Medical Solutions	ActivHeal Dressing	100mm x 100mm	10



Advazorb (T) Plus

- Highly absorbent foam dressing, suitable for those with lots of secretions
- Fenestrated with non-adhesive backing

ltem No.	Supplier	Description	Size	Pack size
CR3900	Insight Medical	Advazorb (T) Plus Dressing	80 mm x 80 mm	20



Advazorb (T) Lite

- Slimline foam dressing, suitable for those with minimal secretions
- Fenestrated with non-adhesive backing

Item No.	Supplier	Description	Size	Pack size
CR4212	Insight Medical	Advazorb (T) Lite Dressing	80 mm x 80 mm	20



Advadraw (T)

- Suitable for low to medium levels of secretions
- Low adherence Reduced trauma on removal

Item No.	Supplier	Description	Size	Pack size
CR4416	Insight Medical	Advadraw (T) Dressing	80 mm x 80 mm	20

Accessories - Swabs and Brushes



Provox© Swab

- Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- The head of the swab is made from foam to avoid loose fibres being left behind
- The white swab easily shows blood or coloured mucus
- Robust and flexible handle for ease of movement inside the tube

ltem No.	Supplier	Description	Size	Pack size
8250	Atos Medical AB	Provox Swab Small	4-5 mm	50
8251	Atos Medical AB	Provox Swab Medium	6-8 mm	50
8252	Atos Medical AB	Provox Swab Large	9-10 mm	50
8258	Atos Medical AB	Provox Swab XtraLarge	11 mm-Above	50



Tracoe[©] Swabs

ltem No.	Supplier	Description	Size	Pack size
TRSWB0150	Atos Medical AB	Trachi-Swab cleaning swab	3.5-5mm	50
TRSWB0250	Atos Medical AB	Trachi-Swab cleaning swab	5-7.5mm	50
TRSWB3050	Atos Medical AB	Trachi-Swab cleaning swab	8mm+	50

Insight Swab

Item No.	Supplier	Description	Size	Pack size
TCS/950	Insight Medical	Insight Swab Small	Up to 6 mm	10
TCS/900	Insight Medical	Insight Swab Large	Up to 9 mm	10

It is recommended that tubes are cleaned with swabs, please check with your healthcare professional for suitability.

Accessories - Wipes and Creams



Provox[©] Cleaning Towel

- A convenient and hygienic way to clean around the stoma
- Small pack size, ideal for care bags, handbags and pockets
- No sting formula
- Does not leave any residue on the skin

Item No.	Supplier	Description	Pack size
7244	Atos Medical AB	Provox Cleaning Towel	10



Provox[©] Skin Barrier

- Leaves a protective film on the skin to help prevent skin irritation
- Use prior to applying adhesive baseplate
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item No.	Supplier	Description	Pack size
8011	Atos Medical AB	Provox Skin Barrier	50



Provox[©] Adhesive Remover

- Helps to remove adhesive comfortably and safely
- Reduces risk of skin irritation
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item No.	Supplier	Description	Pack size
8012	Atos Medical AB	Provox Adhesive Remover	50



Cavilon© Skin Barrier

ltem No.	Supplier	Description	Pack size
3343P	3M Health Care Ltd	Cavilon Skin Barrier Foam Applicator	5
3391G	3M Health Care Ltd	Cavilon Skin Barrier Cream 28g	1
3392G	3M Health Care Ltd	Cavilon Skin Barrier Cream 92g	1
3344E	3M Health Care Ltd	Cavilon Skin Barrier Wipe	30



Medihoney Barrier Cream

ltem No.	Supplier	Description	Pack size
582	Derma Sciences	Medihoney barrier cream 50g	1

Accessories - Shower Protector, Tape, Lubricating Jelly and Irrigation Fluid



Cascade Shower Protector

A flexible, waterproof cover which can be worn over the stoma.

Item No.	Supplier	Description	Pack size
AS3835	Insight Medical	Cascade Shower Protector	2



Item No.SupplierDescriptionSizePack sizeINSIL200Insight MedicalInsil Soft Silicone Perforated Tape20mmx3m1INSIL400Insight MedicalInsil Soft Silicone Perforated Tape40mmx1.5m1



Micropore Tape

Insil Tape

Item No.	Supplier	Description	Size	Pack size
153012	3M Health Care Ltd	Micropore tape	125mm x 5m	1



Optilube™

Item No.	Supplier	Description	Size	Pack size
1120	Optimum Medical Solutions	Optilube Lubricating Jelly Sachets	5g	150 sachets
1121	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	42g	1
1122	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	82g	1



Normasol©

ltem No.	Supplier	Description	Size	Pack size
NOR206B	Molnlycke Health Care	Normasol Sterile Irrigation Fluid	25ml	20



Clinipod©

Irripod©

Item No.	Supplier	Description	Size	Pack size
CLINIPOD	Mayors Healthcare	Clinipod Sodium Chloride 0.9% Irrigation Solution	20ml	25



Item No. Supplier Description Size Pack size IRR24N C D Medical Ltd Irripod Sodium Chloride 0.9% Irrigation Solution 20ml 25

Accessories - Dressing packs and gloves



Sterile Dressing Pack Spec 10

1	Item No.	Supplier	Description	Pack size
N	DP SPEC 10	Generic Supplier	Sterile Dressing Pack Spec 10	12 pcs



Sterile Dressing Pack Spec 35

ltem No.	Supplier	Description	Pack size
DP SPEC 35	Generic Supplier	Sterile Dressing Pack Spec 35	12 pcs



Non-Woven Gauze Swabs

ltem No.	Supplier	Description	Pack size
SWAB SPEC 28	Generic Supplier	Non-Woven Gauze Swabs	100

Gloves



ltem No.	Supplier	Description	Size	Pack size
GLO490B	Generic Supplier	Powder free nitrile gloves	Small	100
GLO489H	Generic Supplier	Powder free nitrile gloves	Medium	100
GLO488L	Generic Supplier	Powder free nitrile gloves	Large	100

Accessories - Protectors



Provox[©] **Protector**

- Thick, absorbent foam stoma protector
- Ideal for those with lots of secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

Item No.	Supplier	Description	Colour	Size	Pack size
7385	Atos Medical AB	Provox Protector	White	Small	10 pcs
7386	Atos Medical AB	Provox Protector	White	Large	10 pcs



Provox[©] **Protector Slim**

- Lightweight, fabric stoma protector
- Easy to breathe through
- Ideal for those with minimal secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

Item No.	Supplier	Description	Colour	Size	Pack size
7387	Atos Medical AB	Provox Protector Slim	White	Small	10 pcs
7388	Atos Medical AB	Provox Protector Slim	Blue	Small	10 pcs
7389	Atos Medical AB	Provox Protector Slim	White	Large	10 pcs
7390	Atos Medical AB	Provox Protector Slim	Blue	Large	10 pcs



Provox[©] **Protector Air**

- Light, foam stoma protector
- Easy to breathe through
- Wipe-able Ideal for those with thick or sticky secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection

ltem No.	Supplier	Description	Colour	Size	Pack size
7391	Atos Medical AB	Provox Protector Air	White	Small	10 pcs
7392	Atos Medical AB	Provox Protector Air	Blue	Small	10 pcs
7393	Atos Medical AB	Provox Protector Air	White	Large	10 pcs
7394	Atos Medical AB	Provox Protector Air	Blue	Large	10 pcs



Sofshield Protector

Item No.	Supplier	Description	Colour	Size	Pack size
AS3698	Insight Medical	Sofshield Protector	White	Small	10 pcs
AS3777	Insight Medical	Sofshield Protector	Blue	Small	10 pcs
AS4056	Insight Medical	Sofshield Protector	Blue	Medium	10 pcs
AS3778	Insight Medical	Sofshield Protector	Beige	Small	10 pcs
AS4063	Insight Medical	Sofshield Protector	Black	Small	10 pcs
AS3669	Insight Medical	Sofshield Protector	White	Large	10 pcs
AS3764	Insight Medical	Sofshield Protector	Blue	Large	10 pcs
AS3765	Insight Medical	Sofshield Protector	Beige	Large	10 pcs
AS4051	Insight Medical	Sofshield Protector	Black	Large	10 pcs



Buchanan Protector

ltem No.	Supplier	Description	Colour	Size	Pack size
BUP0002	Sheffmed Limited	Buchanan Protector	White	Small	10 pcs
BUP0001	Sheffmed Limited	Buchanan Protector	White	Large	10 pcs
LABLT2001	Sheffmed Limited	Buchanan Protector Lite	White	Small	10 pcs
LABLT2002	Sheffmed Limited	Buchanan Protector Lite	Blue	Small	10 pcs
LABLT2003	Sheffmed Limited	Buchanan Protector Lite	Beige	Small	10 pcs
LABLT1001	Sheffmed Limited	Buchanan Protector Lite	White	Large	10 pcs
LABLT1002	Sheffmed Limited	Buchanan Protector Lite	Blue	Large	10 pcs
LABLT1003	Sheffmed Limited	Buchanan Protector Lite	Beige	Large	10 pcs





XtraCare[™] Mini HMEF

XtraCare Mini provides effective filtration along with good humidification for daily protection.

- Effective filtration of bacteria greater than or equal to 99%*
- Effective filtration of viruses greater than or equal to 98%*
- Good humidification
 - Compact design, tailored for children with a trachoestomy
 - Available in three different colors
 - If supplemental oxygen is needed, XtraCare Mini HMEF can be combined with the O2 Adaptor Mini.

ltem No.	Supplier name	Description	Pack size
8004	Atos Medical AB	XtraCare Mini White	30
8005	Atos Medical AB	XtraCare Mini Blue	30
8006	Atos Medical AB	XtraCare Mini Pink	30

*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Freevent® XtraCare Mini HME can never guarantee complete protection. Please read the instructions for use for guidance.

Please check with your healthcare professional for suitability.





O² Adaptor Mini

- Adaptor for administration of supplemental oxygen
- To be used with XtraCare Mini HMEF

Item No.	Supplier name	Description	Pack size
8007	Atos Medical AB	O2 Adaptor Mini	10



TrachPhone® HME

TrachPhone is a versatile and lightweight HME that provides good humidification for different situations.

- A multifunctional HME offering effective lung rehabilitation
- Speaking valve function, press the top lid to occlude the stoma for voicing
- Integrated oxygen supply port (4 mm)
- Integrated suctioning port that allows for tracheal suction without removing the HME

Item No.	Supplier	Description	Pack size
7707	Atos Medical AB	TrachPhone HME	30

Please check with your healthcare professional for suitability.

Kids - Neckbands



Freevent© Neckband One-piece with Twill Ties

ltem No.	Supplier name	Description	Size	Pack size
1651	Atos Medical AB	Freevent Neckband - One-piece	320mm - 420mm	10

Freevent© Neckband Two-piece with VELCRO© fastening



ltem No.	Supplier name	Description	Size	Pack size
1652	Atos Medical AB	Freevent Neckband - Two-piece	230mm - 400mm	10



Trachi-Hold© Tube Holder Small

ltem No.	Supplier	Description	Size	Pack size
TR ACC 0014	Sentient	Trachi-Hold Small	Upto 350mm	10

Trachi-Hold© Tube Holder Mini



ltem No.	Supplier	Description	Size	Pack size
TR ACC 0030	Sentient	Trachi-Hold Mini	Upto 190mm	10

Insight Tube Holder Two-piece



Item No.	Supplier name	Description	Size	Pack size
TH150	Insight Medical	Insight Tube Holder - Two-piece	Upto 345mm	10




Marpac[™] Tube Holder One-piece with Twill Ties

ltem No.	Supplier	Description	Size	Pack size
106	CC Med	Marpac Tube Holder One-Piece	102mm - 178mm	25
107	CC Med	Marpac Tube Holder One-Piece	153mm - 305mm	25

Please check with your healthcare professional for suitability.

$\textbf{Marpac}^{\text{\tiny TM}} \textbf{ Tube Holder Two-piece with VELCRO}^{\texttt{B}} \textbf{/twill hybrid fastener}$





Non-Woven Sterile Gauze Swabs

ltem No.	Supplier name	Description	Pack size
1860A	Generic Supplier	Non-Woven Sterile Gauze Swabs	200



Provox© Swab

- Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- The head of the swab is made from foam to avoid loose fibres being left behind
- The white swab easily shows blood or coloured mucus
- Robust and flexible handle for ease of movement inside the tube

ltem No.	Supplier	Description	Size	Pack size
8250	Atos Medical AB	Provox Swab Small	4-5 mm	50
8251	Atos Medical AB	Provox Swab Medium	6-8 mm	50



Cavilon© Skin Barrier

ltem No.	Supplier	Description	Pack size
3343P	3M Health Care Ltd	Cavilon Skin Barrier Foam Applicator	5
3391G	3M Health Care Ltd	Cavilon Skin Barrier Cream 28g	1
3392G	3M Health Care Ltd	Cavilon Skin Barrier Cream 92g	1
3344E	3M Health Care Ltd	Cavilon Skin Barrier Wipe	30



Optilube™

Item No.	Supplier	Description	Size	Pack size
1120	Optimum Medical Solutions	Optilube Lubricating Jelly Sachets	5g	150 sachets
1121	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	42g	1
1122	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	82g	1

Normasol©



ltem No.	Supplier	Description	Size	Pack size
NOR206B	Molnlycke Health Care	Normasol Sterile Irrigation Fluid	25ml	20



Clinipod©

Item No.	Supplier	Description	Size	Pack size
CLINIPOD	Mayors Healthcare	Clinipod Sodium Chloride 0.9% Irrigation Solution	20ml	25

Irripod©



ltem No.	Supplier	Description	Size	Pack size
IRR24N	C D Medical Ltd	Irripod Sodium Chloride 0.9% Irrigation Solution	20ml	25





Therabite® System Contains:

- 1 x Jaw Mobiliser
- 4 x Bite Pads
- 1 x Hand Aid
- 30 x Range Of Motion Scales
- 1 x Carry Bag
- 1 x Instructions For Use
- 1 x Patient Progress Log

TheraBite®

- For rehabilitation of reduced jaw opening (trismus).
- The TheraBite[®] Jaw Motion Rehabilitation System[™] is a clinically proven, anatomically correct and portable system. Designed to restore mobility after trismus and mandibular hypomobility.
- The TheraBite[®] Jaw Motion Rehabilitation System[™] can also help to reduce muscle pain by reducing joint inflammation. It can prevent reduced muscle strength during radiation treatment (chemotherapy) in the head and neck area.
- The TheraBite Jaw Motion Rehabilitation System is available in both adult and pediatric sizes.

Item No.	Supplier	Description	Quantity
TH001	Atos Medical AB	TheraBite Jaw Motion Reh. System Adult	1
TH002	Atos Medical AB	TheraBite Jaw Motion Reh. System Pediat.	1
PA001	Atos Medical AB	TheraBite Bite Pad, Regular	4
PA002	Atos Medical AB	TheraBite Bite Pad, Edentulous	4
PA003	Atos Medical AB	TheraBite Bite Pad, Pediatric	4

Environmental policies



Fully recyclable packaging

To reduce our carbon footprint we use recyclable materials for every bit of our packaging.



Delivery using electric vehicles

Most of our deliveries are made using electric vehicles which is being expanded over the next 5 years.



Responsible forestry

Where possible, our paper and cardboard are sourced from FSC (Forest Stewardship Council) sources.



Biodegradable materials

Our packaging 'peanuts' are made from biodegradable corn starch and can be dissolved in water.

Contact us

There are several easy ways to get in touch:



Call or text us

Tel: 0800 783 1659 Text: 0753 741 7928



Ē

Email us

info@atos-care.co.uk

Visit our website

www.atos-care.co.uk

Write to us

Atos Medical, Cartwright House, Riverside Business Park, Tottle Road, Nottingham, NG2 1RT

Atos

Atos Care is the distribution arm of Atos Medical, the world-leading supplier of ostomy products for post laryngectomy and tracheostomy patients. Dedicated to making life easier for people living with a neck stoma, Atos Medical has decades of experience, and a deep understanding of what life is like after the operation.

So, although we supply the full range of products, we are much more than just a delivery service. Our expertise and focus mean that we can also provide full professional support and aftercare for you as an individual, wherever you are on your journey. We call this the Atos Circle of Care.

Atos Medical (UK), Cartwright House, Tottle Road, Nottingham, NG2 1RT Tel: 0800 783 1659 Text: 0753 741 7928 Email: info@atos-care.co.uk Web: www.atos-care.co.uk