Care

Laryngectomy Product Catalogue 2025

Bringing the right products and care your way European Contact Centre & Customer Service Awards

2023

SILVER WINNER

Contents

HMEs and Adhesives	12
Provox Life HMEs	12
Provox Life Adhesives	20
Provox Life Tubes	24
Provox Life Buttons	28
Provox	30
Blomsinger System	40
Spirotect System	46

Accessories	50
Tubeholders and Clips	51
Protectors	54
Swabs and Brushes	59
Wipes and Creams	62
ShowerAids	65
Foam Covers	66
Tape and Lubricating Jelly	68
Irrigation Fluid	69
Dressing Packs	70
Gloves	71

Jaw Rehabilitation 72





Welcome to Atos Care

At Atos Care, our mission is to help make life easier for people living with a neck stoma. We understand that having a laryngectomy or tracheostomy can be life changing, and we are here to support you every step of the way.

In this catalogue you will find information on the services available to you and many of the products we supply. If you can't find what you're looking for, please contact our helpful Customer Care team, who can help you find the products you need if they are available on prescription. We supply products made in-house by Atos (the Provox Life, Provox and Tracoe ranges) and products from a variety of other manufacturers, to give you the best choice and selection.

We are here to give you a helping hand, whatever you may need. If you have any questions or feedback, please do not hesitate to get in touch.

How to order

At Atos Care, we are committed to making life easier for people living with a neck stoma. We offer a selection of options for how you can order your supplies; ranging from a monthly service call to a fully automated monthly order.





A friendly monthly service call, email or text

Our helpful Customer Care team can contact you each month to help you place your order and answer any questions you may have. Contact us when you wish to order

If you prefer to contact us when you need to order, you can contact us using our convenient freephone phone line, or by text, email or through our website, if you prefer. Automated convenience order

We can set up a convenient repeat order for you, so that your supplies are automatically ordered and sent to you each month.

Discover our new wesbite **www.atos-care.co.uk** for an easy and convenient way to order products

Order process explained

- 1
- Place an order with our helpful Customer Care team.
- 2

3

- If your order contains new items, we may check suitability with your healthcare professional.
- We contact your GP to request the prescription for your items.*
- 4
- On receiving the prescription**, we send out your supplies directly to your home.

Please allow up to 14 working days from placing your order to receiving your supplies.

The Catalogue doesn't contain everything we can provide, so if there is anything you need but cannot see, then please get in touch with your Customer Care Representative.

^{*}This applies if your GP agrees to receive prescriptions from a third party (i.e. us). If this is not the case, you may be required to place the request directly with your GP surgery.

^{**}Our Electronic Prescriptions Service (EPS) means your GP can send the prescription directly through us through a secure system, so that you don't have to collect it and send it on to us.

Events

We run a range of events, both online and in-person. Our MyLife events are one of the ways we can help support you to breathe better, whatever you do.

We hold a range of helpful webinars, practical workshops and interactive experiences. Our MyLife events include helpful presentations about relevant topics, as well as a chance to ask questions and hear from other people who know and understand about laryngectomy.

Access to experts

Supported events can also include healthcare and other professionals who can share valuable information, tips and best practice that might be hard to find on your own. Past events have included presentations and speakers on themes such as:


- Looking after your lungs
- Travel
- Breathing exercises
- Dealing with the unexpected



Atos MyLife app

Did you know that we have a helpful App, full of informative videos, tutorials and tips about living with a laryngectomy?

Featuring real patient stories as well as presentations from our in-house Speech Therapists and guest experts.



The App is divided into three main sections making it easy for you to find what you are looking for:

Inspiration

Articles and original videos with patients and clinicians sharing advice and unique stories. Topics include breathing, speaking, care and daily life.

Exercises

A tool to help you practice voice and breathing exercises through a series of videos offering useful tips to help you overcome any challenges.

Products

A presentation of the Provox Life products with information to help understand how to use them best within your daily life.



Award winning app

We are thrilled to announce that our MyLife app won four awards at the 12th annual Lovie awards.

- Gold Lovie Winner: Utilities
 and Services
- Silver Lovie Winner: Wellness & Lifestyle
- Bronze Lovie Winner: Best
 practices
 - People's Lovie Winner: Best practices

Download the **MyLife app** now











Provox[®] Life[™]



Provox[®] Life[™] Home



This HME is designed to improve your lung health when you are taking it easy at home. It features optimised humidification and breathability for healthier lungs.

- Excellent humidification
- Comfortable breathability
- Low profile and discreet look

- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis



Provox[®] Life[™] Go

This HME provides lower breathing resistance, helping make it easier to breathe when you are out and about and participating in daily activities.

- Good breathability
- Good humidification
- Low-profile and discreet look
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis

ltem No.	Supplier	Description	Pack size
8311	Atos Medical AB	Provox Life Home HME	30 рс
8310	Atos Medical AB	Provox Life Go HME	30 рс



Provox® Life™ Night

This HME helps reduce coughing for a more restful nights sleep. It offers our highest humidification which helps reduce mucus production and coughing, allowing you to rest whilst improving lung health.

- Excellent humidification
- Comfortable breathability for nighttime breathing
- Reduces mucus production and coughing
- Soft and smooth for nighttime comfort
- Side openings to prevent occlusion while sleeping



Provox[®] Life[™] Energy

This HME has good humidification and low breathing resistance, which is perfect for when you are physically active.

- High breathability
- Good humidification
- Sporty look
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis
- Front vents to prevent accidental occlusion during deep breathing

ltem	No.	Supplier	Description	Pack size
8262	2	Atos Medical AB	Provox Life Night HME	30 pc
8312		Atos Medical AB	Provox Life Energy HME	30 рс



Provox® Life™ Protect



This HME has a two-way electrostatic filtration that effectively filters bacteria, viruses, dust, and pollen. It also has a comfortable breathing resistance that warms and humidifies the air you breathe.

- A highly effective bacterial and viral filter (greater than or equal to 98%)*
- Protects you from viruses, bacteria, dust and pollen
- Reduces mucus production and coughing
- Easy to press for speaking with a voice prosthesis

Item No.	Supplier	Description	Pack size
8313	Atos Medical AB	Provox Life Protect HME	15 pc

Ready for an upgrade? Text 'LIFE' to 07537 417 928 or call 0800 783 1659

*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Provox Life™ Protect HME can never guarantee complete protection. Please read the instructions for use for guidance.



Provox® Life™ FreeHands



Provox[®] FreeHands HME offers both the freedom to speak hands-free and with finger occlusion. For use with a voice prosthesis. To be used with Provox[®] FreeHands FlexiVoice[™].

- Reduces mucus production and coughing
- Easy removal of HME without the need for a removal aid
- Why not try with the Provox Life[™] Stability Adhesive for extra confidence when speaking with Provox Life[™] FreeHands

ltem No.	Supplier	Description	Pack size
7440	Atos Medical AB	Provox Life FreeHands HME	30 рс

"I am very happy with the Provox Life™ Stability adhesive and don't feel nervous it's going to blow off when voicing."

- Provox Life[™] User



Provox® FreeHands FlexiVoice™

Gives the ability to speak hands-free. For use with a voice prosthesis.

- Speak hands-free for more natural communication
- Can be covered with a finger if needed
- Quick and easy to set up
- Helpful arch accessory helps keep clothing away

- Available in four strengths to suit your needs
- Compatible with Provox Life[™] Adhesives using Provox Life[™] FreeHands HME, or with Provox[®] Adhesives using Provox[®] FreeHands HMEs

Item No.	Supplier	Description	Pack size
8161	Atos Medical AB	Provox FreeHands FlexiVoice Light	1 рс
8162	Atos Medical AB	Provox FreeHands FlexiVoice Medium	1 рс
8163	Atos Medical AB	Provox FreeHands FlexiVoice Strong	1 рс
8166	Atos Medical AB	Provox FreeHands FlexiVoice XtraStrong	1 рс

Your clinician may wish to do a pressure test with you to determine the best FlexiVoice™ for you to use.



Provox[®] FreeHands Support[™]

To help reduce stoma movement, improve voice quality and prevent loss of baseplate seal while speaking hands-free.

- Helps keep a good seal with your baseplate while speaking hands-free
- Three models available to suit different stoma depths
- Adhesive panel fixes to the chest to reduce baseplate movement whilst speaking

Item No.	Supplier	Description	Pack size
8021	Atos Medical AB	Provox FreeHands Support Flat	1 рс
8022	Atos Medical AB	Provox FreeHands Support Medium	1 рс
8023	Atos Medical AB	Provox FreeHands Support Deep	1 рс
8024	Atos Medical AB	Provox FreeHands Support Adhesive	15 pcs

Provox® Life™ Standard Adhesive

An everyday adhesive that provides a strong seal around your stoma. A soft, flexible and low-profile adhesive that forms a strong seal and provides a secure fit.

- Strong seal
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit[™] technology comfortable low-profile connector
- Compatible with all Provox Life[™] HMEs
- Ideal for medium to flat stomas

ltem No.	Supplier	Description	Pack size
7460	Atos Medical AB	Provox Life Standard Adhesive - Round	30 рс
7461	Atos Medical AB	Provox Life Standard Adhesive - Oval	30 рс
7462	Atos Medical AB	Provox Life Standard Adhesive - Plus	30 рс



\ Provox® Life™ Sensitive Adhesive

A gentle everyday adhesive for sensitive skin. A comfortable, gentle and low-profile adhesive that forms a good seal for a secure fit.

- Good seal
- Hypoallergenic material for sensitive skin
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit[™] technology comfortable low-profile connector
- Ideal for sensitive skin
- Compatible with All Provox Life[™] HMEs

ltem No.	Supplier	Description	Pack size
7463	Atos Medical AB	Provox Life Sensitive Adhesive - Round	30 рс
7464	Atos Medical AB	Provox Life Sensitive Adhesive - Oval	30 рс
7466	Atos Medical AB	Provox Life Sensitive Adhesive - Plus	30 рс









∧ Provox® Life™ Stability Adhesive

0

A firm everyday adhesive that provides a secure attachment, especially if your stoma is deep. It is ideal for those who speak frequently or speak hands-free.

- Strong seal
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit[™] technology comfortable low-profile connector
- Compatible with all Provox Life[™] HMEs

- Ideal for deep stomas or speaking handsfree
- Reinforced centre panel for extra stability

Item No.	Supplier	Description	Pack size
8263	Atos Medical AB	Provox Life Stability Adhesive	15 pc

∧ Provox® Life™ Night Adhesive

This adhesive allows your skin to recover while you sleep or during the day. A skin-friendly hydrogel adhesive that soothes and cools the skin.

- Soft and smooth for nighttime comfort
- Hydrogel material helps soothe and cool the skin while you sleep
- Innovative clover-shaped adhesive for easy application and better fit
- SecureFit[™] technology comfortable low-profile connector
- Compatible with all Provox Life[™] HMEs
- Ideal for sensitive, irritated skin
- Can be worn during the day when needed

ltem No.	Supplier	Description	Pack size
8261	Atos Medical AB	Provox Life Night Adhesive	15 рс



Provox[©] Life[™] LaryTube Standard

A flexible silicone tube to attach a Provox Life[™] HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox© TubeHolder or Provox© LaryClips.

Item No.	Supplier	Description	ID	OD	L	Pack size
7409	Atos Medical AB	Provox Life LaryTube 8/27 - Standard	9.5	12	27	1 рс
7410	Atos Medical AB	Provox Life LaryTube 8/36 - Standard	9.5	12	36	1 рс
7411	Atos Medical AB	Provox Life LaryTube 8/55 - Standard	9.5	12	55	1 рс
7412	Atos Medical AB	Provox Life LaryTube 9/27 - Standard	10.5	13.5	27	1 рс
7413	Atos Medical AB	Provox Life LaryTube 9/36 - Standard	10.5	13.5	36	1 рс
7414	Atos Medical AB	Provox Life LaryTube 9/55 - Standard	10.5	13.5	55	1 рс
7415	Atos Medical AB	Provox Life LaryTube 10/27 - Standard	12	15	27	1 рс
7416	Atos Medical AB	Provox Life LaryTube 10/36 - Standard	12	15	36	1 рс
7417	Atos Medical AB	Provox Life LaryTube 10/55 - Standard	12	15	55	1 рс
7418	Atos Medical AB	Provox Life LaryTube 12/27 - Standard	13.5	17	27	1 рс
7419	Atos Medical AB	Provox Life LaryTube 12/36 - Standard	13.5	17	36	1 рс
7420	Atos Medical AB	Provox Life LaryTube 12/55 - Standard	13.5	17	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm)

L = Length (mm)

Please check with your healthcare professional for suitability.



Provox[©] Life[™] LaryTube Standard with Ring

A flexible silicone tube to attach a Provox Life[™] HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with $\mathsf{Provox}^{\mathbb{C}}$ Life^t Adhesives.

ltem No.	Supplier	Description	ID	OD	L	Pack size
7421	Atos Medical AB	Provox Life LaryTube 8/36 - Std with Ring	9.5	12	36	1 рс
7422	Atos Medical AB	Provox Life LaryTube 8/55 - Std with Ring	9.5	12	55	1 рс
7423	Atos Medical AB	Provox Life LaryTube 9/36 - Std with Ring	10.5	13.5	36	1 рс
7424	Atos Medical AB	Provox Life LaryTube 9/55 - Std with Ring	10.5	13.5	55	1 рс
7425	Atos Medical AB	Provox Life LaryTube10/36 - Std with Ring	12	15	36	1 рс
7426	Atos Medical AB	Provox Life LaryTube 10/55 - Std with Ring	12	15	55	1 рс
7427	Atos Medical AB	Provox Life LaryTube 12/36 - Std with Ring	13.5	17	36	1 рс
7428	Atos Medical AB	Provox Life LaryTube 12/55 - Std with Ring	13.5	17	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)



Provox[©] Life[™] LaryTube Fenestrated

A flexible silicone tube with holes to attach a Provox Life[™] HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox[©] TubeHolder or Provox[©] LaryClips.

Item No.	Supplier	Description	ID	OD	L	Pack size
7429	Atos Medical AB	Provox Life LaryTube 8/36 - Fenestrated	9.5	12	36	1 рс
7430	Atos Medical AB	Provox Life LaryTube 8/55 - Fenestrated	9.5	12	55	1 рс
7431	Atos Medical AB	Provox Life LaryTube 9/36 - Fenestrated	10.5	13.5	36	1 рс
7432	Atos Medical AB	Provox Life LaryTube 9/55 - Fenestrated	10.5	13.5	55	1 рс
7433	Atos Medical AB	Provox Life LaryTube 10/36 - Fenestrated	12	15	36	1 рс
7434	Atos Medical AB	Provox Life LaryTube 10/55 - Fenestrated	12	15	55	1 рс
7435	Atos Medical AB	Provox Life LaryTube 12/36 - Fenestrated	13.5	17	36	1 рс
7436	Atos Medical AB	Provox Life LaryTube 12/55 - Fenestrated	13.5	17	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)



Provox[©] Life[™] LaryTube Fenestrated with Ring

A flexible silicone tube with holes to attach a Provox Life[™] HME or accessory, while making sure the stoma maintains its size and shape.

Can be used with or without a voice prosthesis. To be secured with Provox© Life™ Adhesives.

ltem No.	Supplier	Description	ID	OD	L	Pack size
8048	Atos Medical AB	Provox Life LaryTube 8/36 - Fen with Ring	9.5	12	36	1 рс
8049	Atos Medical AB	Provox Life LaryTube 8/55 - Fen with Ring	9.5	12	55	1 рс
8050	Atos Medical AB	Provox Life LaryTube 9/36 - Fen with Ring	10.5	13.5	36	1 рс
8051	Atos Medical AB	Provox Life LaryTube 9/55 - Fen with Ring	10.5	13.5	55	1 рс
8052	Atos Medical AB	Provox Life LaryTube 10/36 - Fen with Ring	12	15	36	1 рс
8053	Atos Medical AB	Provox Life LaryTube 10/55 - Fen with Ring	12	15	55	1 рс
8054	Atos Medical AB	Provox Life LaryTube 12/36 - Fen with Ring	13.5	17	36	1 рс
8055	Atos Medical AB	Provox Life LaryTube 12/55 - Fen with Ring	13.5	17	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm)

L = Length (mm)



Provox[©] Life[™] LaryButton

A soft silicone stoma button to attach a Provox Life $^{\rm \tiny M}$ HME or accessory, while keeping the stoma stable in size.

To be secured with Provox[©] Tube holder or Provox[©] LaryClips.

ltem No.	Supplier	Description	OD	L	Pack size
8040	Atos Medical AB	Provox Life LaryButton 12/8	12	8	1 рс
8041	Atos Medical AB	Provox Life LaryButton 12/18	12	18	1 рс
8042	Atos Medical AB	Provox Life LaryButton 14/8	14	8	1 рс
8043	Atos Medical AB	Provox Life LaryButton 14/18	14	18	1 рс
8044	Atos Medical AB	Provox Life LaryButton 16/8	16	8	1 рс
8045	Atos Medical AB	Provox Life LaryButton 16/18	16	18	1 рс
8046	Atos Medical AB	Provox Life LaryButton 18/8	18	8	1 рс
8047	Atos Medical AB	Provox Life LaryButton 18/18	18	18	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm) "I have actually been getting decent sleep since the first night."

- Provox Life™ User





Provox[©] XtraMoist[™] HME

For excellent humidification when relaxing to help reduce coughing and mucus.

Item No.	Supplier	Description	Pack size
7290	Atos Medical AB	Provox XtraMoist HME	30 pcs



Provox[©] XtraFlow[™] HME

For humidification and lower breathing resistance when being active.

Item No.	Supplier	Description	Pack size
7291	Atos Medical AB	Provox XtraFlow HME	30 pcs



Provox© Micron[™] HME

A highly effective virus and bacteria filter (over 99%*) with good humidification.

ltem No.	Supplier	Description	Pack size
7248	Atos Medical AB	Provox Micron HME	30 pcs



Provox© Luna™ HME

A soft HME for excellent humidification when sleeping.

This HME should be used with Provox Luna adhesive.

ltem No.	Supplier	Description	Pack size
8013	Atos Medical AB	Provox Luna HME	30 pcs



Provox[®] FreeHands HME[®] Moist

Use with Provox FreeHands FlexiVoice for humidification when relaxing.

ltem No.	Supplier	Description	Pack size
8220	Atos Medical AB	Provox FreeHands HME Moist	30 pcs



Provox® FreeHands HME® Flow

Use with Provox FreeHands FlexiVoice for a lower breathing resistance when being active.

ltem No.	Supplier	Description	Pack size
8221	Atos Medical AB	Provox FreeHands HME Flow	30 pcs

*Please note: Since pathogens can enter and leave the human body in other ways (such as the mouth, nose, and eyes), Provox Micron HME can never guarantee complete protection. Please read the instructions for use for guidance.



Provox[©] **FlexiDerm[™] Adhesive**

Gives you a good seal around your stoma.

Item No.	Supplier	Description	Pack size
7253	Atos Medical AB	Provox FlexiDerm Adhesive - Round	20 pcs
7254	Atos Medical AB	Provox FlexiDerm Adhesive - Oval	20 pcs
7331	Atos Medical AB	Provox FlexiDerm Adhesive - Plus	20 pcs



Provox[©] **OptiDerm**[™] **Adhesive**

An adhesive if your skin feels sore or irritated.

ltem No.	Supplier	Description	Pack size
7255	6 Atos Medical AB Provox OptiDerm Adhesive - Oval		20 pcs
7256			20 pcs
7332			20 pcs



Provox[©] StabiliBase[™] Adhesive

Gives a superior seal and support when speaking hands-free.

Acrylic adhesive with firm concave base, providing support and strong adhesion.

ltem No.	Supplier	Description	Pack size
7289	Atos Medical AB	Provox StabiliBase Adhesive	15 pcs



Provox[©] **StabiliBase OptiDerm[™] Adhesive**

When you have a deep stoma and sensitive skin.

ltem No.	Supplier	Description	Pack size
7318	Atos Medical AB	Provox StabiliBase OptiDerm Adhesive	15 pcs



Provox© XtraBase™ Adhesive

For speaking with finger occlusion or hands-free.

Strong, acrylic adhesive with reinforced concave base for extra support.

ltem No.	Supplier	Description	Pack size
7265	Atos Medical AB	Provox XtraBase Adhesive	20 pcs



Provox© Luna™ Adhesive

Designed to help you sleep comfortably and soothe your skin.

ltem No.	Supplier	Description	Pack size
8014	Atos Medical AB	Provox Luna Adhesive	15 pcs



Provox[©] Adhesive Strip

Protects Provox Life[™] Night Adhesive or Provox[®] Luna Adhesive when showering.

Item No.	Supplier	Description	Pack size
8015	Atos Medical AB	Provox Adhesive Strip	10 pcs



Provox[©] LaryTube Standard

A soft silicone tube used to keep the stoma stable in size.

Can be used with or without a voice prosthesis. To be secured with Provox© TubeHolder or Provox© LaryClips.

ltem No.	Supplier	Description	ID	OD	L	Pack size
7601	Atos Medical AB	Provox LaryTube 8/27 - Standard	9.5	12	27	1 рс
7602	Atos Medical AB	Provox LaryTube 8/36 - Standard	9.5	12	36	1 рс
7603	Atos Medical AB	Provox LaryTube 8/55 - Standard	9.5	12	55	1 рс
7605	Atos Medical AB	Provox LaryTube 9/27 - Standard	10.5	13.5	27	1 рс
7606	Atos Medical AB	Provox LaryTube 9/36 - Standard	10.5	13.5	36	1 рс
7607	Atos Medical AB	Provox LaryTube 9/55 - Standard	10.5	13.5	55	1 рс
7609	Atos Medical AB	Provox LaryTube 10/27 - Standard	12	15	27	1 рс
7610	Atos Medical AB	Provox LaryTube 10/36 - Standard	12	15	36	1 рс
7611	Atos Medical AB	Provox LaryTube 10/55 - Standard	12	15	55	1 рс
7613	Atos Medical AB	Provox LaryTube 12/27 - Standard	13.5	17	27	1 рс
7614	Atos Medical AB	Provox LaryTube 12/36 - Standard	13.5	17	36	1 рс
7615	Atos Medical AB	Provox LaryTube 12/55 - Standard	13.5	17	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)

Please check with your healthcare professional for suitability.



Provox[©] LaryTube Standard with Ring

A soft silicone tube used to keep the stoma stable in size.

Can be used with or without a voice prosthesis. Attaches directly into Provox© adhesive baseplates; no need for neckbands.

ltem No.	Supplier	Description	ID	OD	L	Pack size
7624	Atos Medical AB	Provox LaryTube 8/36 - Std with Ring	9.5	12	36	1 рс
7625	Atos Medical AB	Provox LaryTube 8/55 - Std with Ring	9.5	12	55	1 рс
7626	Atos Medical AB	Provox LaryTube 9/36 - Std with Ring	10.5	13.5	36	1 рс
7627	Atos Medical AB	Provox LaryTube 9/55 - Std with Ring	10.5	13.5	55	1 рс
7628	Atos Medical AB	Provox LaryTube10/36 - Std with Ring	12	15	36	1 рс
7629	Atos Medical AB	Provox LaryTube 10/55 - Std with Ring	12	15	55	1 рс
7630	Atos Medical AB	Provox LaryTube 12/36 - Std with Ring	13.5	17	36	1 рс
7631	Atos Medical AB	Provox LaryTube 12/55 - Std with Ring	13.5	17	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)



Provox[©] LaryTube Fenestrated

A soft silicone tube used to keep the stoma stable in size.

Small holes (fenestrations) are designed for use with a voice prosthesis. To be secured with Provox[©] TubeHolder or Provox[©] LaryClips.

ltem No.	Supplier	Description	ID	OD	L	Pack size
7637	Atos Medical AB	Provox LaryTube 8/36 - Fenestrated	9.5	12	36	1 pc
7638	Atos Medical AB	Provox LaryTube 8/55 - Fenestrated	9.5	12	55	1 pc
7640	Atos Medical AB	Provox LaryTube 9/36 - Fenestrated	10.5	13.5	36	1 pc
7641	Atos Medical AB	Provox LaryTube 9/55 - Fenestrated	10.5	13.5	55	1 рс
7643	Atos Medical AB	Provox LaryTube 10/36 - Fenestrated	12	15	36	1 pc
7644	Atos Medical AB	Provox LaryTube 10/55 - Fenestrated	12	15	55	1 рс
7646	Atos Medical AB	Provox LaryTube 12/36 - Fenestrated	13.5	17	36	1 рс
7647	Atos Medical AB	Provox LaryTube 12/55 - Fenestrated	13.5	17	55	1 pc

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)



Provox[©] LaryButton

A soft silicone stoma button used to keep the stoma stable in size. To be secured with Provox[©] TubeHolder or Provox[©] LaryClips.

Item No.	Supplier	Description	OD	L	Pack size
7671	Atos Medical AB	Provox LaryButton 12/8	12	8	1 рс
7672	Atos Medical AB	Provox LaryButton 14/8	14	8	1 рс
7673	Atos Medical AB	Provox LaryButton 16/8	16	8	1 рс
7674	Atos Medical AB	Provox LaryButton 18/8	18	8	1 рс
7685	Atos Medical AB	Provox LaryButton 12/18	12	18	1 рс
7686	Atos Medical AB	Provox LaryButton 14/18	14	18	1 рс
7687	Atos Medical AB	Provox LaryButton 16/18	16	18	1 рс
7688	Atos Medical AB	Provox LaryButton 18/18	18	18	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)

Please check with your healthcare professional for suitability.

Blom-Singer©



Blom-Singer© SpeakFree© HME

Single-use adjustable hands-free valve.

Item No.	Supplier	Description	Pack size
BE 1090 EF	Severn Healthcare Technologies	Blom-Singer SpeakFree HME EasyFlow	30 pcs
BE 1090 EZ	Severn Healthcare Technologies	Blom-Singer SpeakFree HME ClassicFlow	30 pcs



Blom-Singer© Day and Night© ClassicFlow HME

Touch to speak HME.

ltem No.	Supplier	Description	Pack size
BE 1091	Severn Healthcare Technologies	Blom-Singer Day and Night ClassicFlow	30 pcs



Blom-Singer© Day and Night© EasyFlow HME

Touch to speak HME.

ltem No.	Supplier	Description	Pack size
BE 1092	Severn Healthcare Technologies	Blom-Singer Day and Night EasyFlow	30 pcs



Blom-Singer© AccuFit© Adhesive

Flexible, translucent adhesive.

Item No.	Supplier	Description	Pack size
BE 6082	Severn Healthcare Technologies	Blom-Singer AccuFit Adhesive - Round	30 pcs
BE 6083	Severn Healthcare Technologies	Blom-Singer AccuFit Adhesive - Oval	30 pcs
BE 6084	Severn Healthcare Technologies	Blom-Singer AccuFit Adhesive - Oval Extra	30 pcs



Blom-Singer© HydroFit© Adhesive

Hydrocolloid adhesive suitable for sensitive skin.

ltem No.	Supplier	Description	Pack size
BE 6085	Severn Healthcare Technologies	Blom-Singer HydroFit Hydrocolloid Adhesive - Round	30 pcs
BE 6086	Severn Healthcare Technologies	Blom-Singer HydroFit Hydrocolloid Adhesive - Oval	30 pcs
BE 6087	Severn Healthcare Technologies	Blom-Singer HydroFit Hydrocolloid Adhesive - Oval Extra	30 pcs

Please check with your healthcare professional for suitability.



Blom-Singer© Tracheostoma Valve Housing

A reusable silicone housing for use with Tape Adhesive discs. Fits with the ATSVII or Humidfilter Holder.



ltem No.	Supplier	Description	Pack size
BE6038	Severn Healthcare Technologies	Blom-Singer Tracheostoma Valve Housing - Standard	1 рс
BE6039	Severn Healthcare Technologies	Blom-Singer Tracheostoma Valve Housing - Large	1 рс



Blom-Singer© Tape Adhesive Discs

Double-sided adhesive tape discs for use with the Tracheostoma Valve Housing.

Item No.	Supplier	Description	Pack size
BE 6041	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Standard, Regular	30 pcs
BE 6042	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Large, Regular	30 pcs
BE 6043	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Standard, Thin	30 pcs
BE 6044	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Large, Thin	30 pcs
BE6035	Severn Healthcare Technologies	Blom-Singer Tape Adhesive Discs - Large, Heavy Duty	30 pcs



Blom-Singer© Stomasoft© Tube

Soft and flexible laryngectomy tube.

Item No.	Supplier	Description	ID	OD	L	Pack size
BE 6698	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Standard	9.5	12	36	1 рс
BE 6699	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Standard	9.5	12	55	1 рс
BE 6700	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Standard	10.5	13.5	36	1 рс
BE 6701	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Standard	10.5	13.5	55	1 рс
BE 6702	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Standard	12	15	36	1 рс
BE 6703	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Standard	12	15	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)



Blom-Singer© Stomasoft© Tube

Soft and flexible laryngectomy tube with fenestration.

Item No.	Supplier	Description	ID	OD	L	Pack size
BE 6698F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Fenestrated	9.5	12	36	1 рс
BE 6699F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 8 - Fenestrated	9.5	12	55	1 рс
BE 6700F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Fenestrated	10.5	13.5	36	1 рс
BE 6701F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 9 - Fenestrated	10.5	13.5	55	1 рс
BE 6702F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Fenestrated	12	15	36	1 рс
BE 6703F	Severn Healthcare Technologies	Blom-Singer Stomasoft Tube Size 10 - Fenestrated	12	15	55	1 рс

ID = Inner diameter (mm) OD = Outer diameter (mm) L = Length (mm)

SpiroTect©



SpiroTect® Touch2Speak HME

Easily obtain airflow closure and speech by using a light finger touch.

Item No.	Supplier name	Description	Pack size
PTC1004	Sheffmed Limited	Spirotect Touch2Speak HME	30 pcs



SpiroTect® Push2Speak HME

Easily obtain airflow closure and speech by gently pushing the top of the cassette.

ltem No.	Supplier name	Description	Pack size
PTC1005	Sheffmed Limited	Spirotect Push2Speak HME	30 pcs



SpiroTect[®] Fabric Adhesive

Fabric dressings that secure HME cassettes.

ltem No.	Supplier name	Description	Pack size
PTW100	Sheffmed Limited	Spirotect Fabric Adhesive - Oval	20 pcs
PTW101	Sheffmed Limited	Spirotect Fabric Adhesive - Round	20 pcs



SpiroTect[®] Sensitive Adhesive

Made from hydrocolloid that secure HME cassettes.

ltem No.	Supplier name	Description	Pack size
PTW200	Sheffmed Limited	Spirotect Sensitive Adhesive - Oval	20 pcs
PTW201	Sheffmed Limited	Spirotect Sensitive Adhesive - Round	20 pcs

0

SpiroTect[®] Clear & Strong Adhesive

Secure HME cassettes at the stoma site.

ltem	n No.	Supplier name	Description	Pack size
PTV	V300	Sheffmed Limited	Spirotect Clear & Strong Adhesive - Oval	20 pcs
PTV	V301	Sheffmed Limited	Spirotect Clear & Strong Adhesive - Round	20 pcs

SpiroTect[®] Extra Sticky Adhesive

Developed with extra adhesion to secure HME cassettes.

ltem No.	Supplier name	Description	Pack size
PTW501	Sheffmed Limited	Spirotect Extra Sticky Adhesive - Round	20 pcs



SpiroTect[®] Secure Adhesive

Supportive and secure positioning around the stoma site.

ltem No.	Supplier name	Description	Pack size
PTW700	Sheffmed Limited	Spirotect Secure Adhesive - Convex	20 pcs

Accessories

50

Accessories - TubeHolders and Clips



Provox© TubeHolder

This can be used to hold a laryngectomy tube or button in place. The discreet clips fit into the connectors of a Provox[®] or Provox Life[™] LaryTube and both Provox[®] and Provox Life[™] LaryButtons.

- Discreet colour
- Side closure with velcro for adjustable length
- High quality material with slight stretch for comfort
- Compatible with Provox Life™

Item No.	Supplier	Description	Size	Pack size
7668	Atos Medical AB	Provox TubeHolder	Adjustable	1



Can be used to hold a laryngectomy tube or button in place. The discreet clips fit into the connectors of a Provox[®] and Provox Life[™] LaryTube or LaryButton.

- Discreet colour
- Adhesive backing attaches securely to skin
- Easy fastenings for flexible positioning
- Compatible with Provox Life™

Item No.	Supplier	Description	Size	Pack size
7669	Atos Medical AB	Provox LaryClip	Adjustable	20 pads/8 clips



Freevent® Neckband One-piece with Twill Ties

- · One-piece neckband with twill ties, designed for extra security
- Simply cut and adjust to the exact size needed
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort

Item No.	Supplier	Description	Size	Pack size
1661	Atos Medical AB	Freevent Neckband - One-piece	Upto 550mm	10



Freevent® Neckband Two-piece with VELCRO® fastening

- High quality adjustable neckband with VELCRO® fastening designed to provide secure tube positioning
- Simple application with VELCRO® fastening, adjustable to different neck sizes
- Made of soft foam to help minimise skin irritation
- Rounded edges for comfort

Item No.	Supplier	Description	Size	Pack size
1662	Atos Medical AB	Freevent Neckband - Two-piece	Upto 550mm	10



Insight Tube Holder

ltem No.	Supplier	Description	Size	Pack size
TH100	Insight Medical	Insight Tube Holder - Two-piece	Upto 470mm	10

Accessories - Protectors



Provox[©] **Protector**

- Thick, absorbent foam stoma protector
- Ideal for those with lots of secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

ltem No.	Supplier	Description	Colour	Size	Pack size
7385	Atos Medical AB	Provox Protector	White	Small	10 pcs
7386	Atos Medical AB	Provox Protector	White	Large	10 pcs



Provox[©] **Protector Slim**

- Lightweight, fabric stoma protector
- Easy to breathe through
- Ideal for those with minimal secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection, especially in winter

ltem No.	Supplier	Description	Colour	Size	Pack size
7387	Atos Medical AB	Provox Protector Slim	White	Small	10 pcs
7388	Atos Medical AB	Provox Protector Slim	Blue	Small	10 pcs
7389	Atos Medical AB	Provox Protector Slim	White	Large	10 pcs
7390	Atos Medical AB	Provox Protector Slim	Blue	Large	10 pcs



Provox[©] **Protector Air**

- Light, foam stoma protector
- Easy to breathe through
- Wipe-able Ideal for those with thick or sticky secretions
- Comfortable and discreet
- Soft-lock mechanism for comfort
- Hand washable up to 3 times
- Suitable to be worn over an HME for additional protection

Item No.	Supplier	Description	Colour	Size	Pack size
7391	Atos Medical AB	Provox Protector Air	White	Small	10 pcs
7392	Atos Medical AB	Provox Protector Air	Blue	Small	10 pcs
7393	Atos Medical AB	Provox Protector Air	White	Large	10 pcs
7394	Atos Medical AB	Provox Protector Air	Blue	Large	10 pcs



Sofshield Protector

ltem No.	Supplier	Description	Colour	Size	Pack size
AS3698	Insight Medical	Sofshield Protector	White	Small	10 pcs
AS3777	Insight Medical	Sofshield Protector	Blue	Small	10 pcs
AS4056	Insight Medical	Sofshield Protector	Blue	Medium	10 pcs
AS3778	Insight Medical	Sofshield Protector	Beige	Small	10 pcs
AS4063	Insight Medical	Sofshield Protector	Black	Small	10 pcs
AS3669	Insight Medical	Sofshield Protector	White	Large	10 pcs
AS3764	Insight Medical	Sofshield Protector	Blue	Large	10 pcs
AS3765	Insight Medical	Sofshield Protector	Beige	Large	10 pcs
AS4051	Insight Medical	Sofshield Protector	Black	Large	10 pcs



Buchanan Protector

ltem No.	Supplier	Description	Colour	Size	Pack size
BUP0002	Sheffmed Limited	Buchanan Protector	White	Small	10 pcs
BUP0001	Sheffmed Limited	Buchanan Protector	White	Large	10 pcs
LABLT2001	Sheffmed Limited	Buchanan Protector Lite	White	Small	10 pcs
LABLT2002	Sheffmed Limited	Buchanan Protector Lite	Blue	Small	10 pcs
LABLT2003	Sheffmed Limited	Buchanan Protector Lite	Beige	Small	10 pcs
LABLT1001	Sheffmed Limited	Buchanan Protector Lite	White	Large	10 pcs
LABLT1002	Sheffmed Limited	Buchanan Protector Lite	Blue	Large	10 pcs
LABLT1003	Sheffmed Limited	Buchanan Protector Lite	Beige	Large	10 pcs

Accessories - Swabs and Brushes



Provox© Swab

- Strong, sturdy, smooth swab head, designed to clean the surface without scratching
- The head of the swab is made from foam to avoid loose fibres being left behind
- The white swab easily shows blood or coloured mucus
- Robust and flexible handle for ease of movement inside the tube

Item No.	Supplier	Description	Size	Pack size
8250	Atos Medical AB	Provox Swab Small	4-5 mm	50
8251	Atos Medical AB	Provox Swab Medium	6-8 mm	50
8252	Atos Medical AB	Provox Swab Large	9-10 mm	50
8258	Atos Medical AB	Provox Swab XtraLarge	11 mm-Above	50



Insight Swab

Item No.	Supplier	Description	Size	Pack size
TCS/950	Insight Medical	Insight Swab Small	Up to 6 mm	10
TCS/900	Insight Medical	Insight Swab Large	Up to 9 mm	10



Provox© TubeBrush

- Easy to grip handle
- Bristles that extend further down the wire of the brush
- Softer bristles and a soft cotton tip to help reduce the risk of damage to the surface of the tube

ltem No.	Supplier	Description	Size	Pack size
7661	Atos Medical AB	Provox TubeBrush	12 mm	6
7660	Atos Medical AB	Provox TubeBrush	8 mm	6



Provox[©] Voice Prosthesis Brush

- Rounded tip to protect against accidental injury
- Safety wings to prevent unlikely inhalation of the Provox Brush
- Short handle for easier handling
- Use the Provox Brush to insert the Provox Plug

ltem No.	Supplier	Description	Size	Pack size
7204	Atos Medical AB	Provox Brush	4.5-10 mm	6
7225	Atos Medical AB	Provox Brush XL	12-18 mm	6
8404	Atos Medical AB	Provox Brush Long	4-10 mm	6
8425	Atos Medical AB	Provox Brush XL Long	12-15 mm	6

Please check with your healthcare professional for suitability.



Blom-Singer© Brush

Item No.	Supplier	Description	Size	Pack size
IN4064	Severn Healthcare Technologies	Blom-Singer Brush	4-8 mm	3
IN4065	Severn Healthcare Technologies	Blom-Singer Brush	10-16 mm	3
IN4066	Severn Healthcare Technologies	Blom-Singer Brush	18 mm-Above	3

Insight Brush

Item No.	Supplier	Description	Size	Pack size
PB/50	Insight Medical	Insight Brush	One size	10



Kapitex Voiceline Brush

ltem No.	Supplier	Description	Size	Pack size
LAPCB 0003	Kapitex Healthcare	Voiceline brush	One size	10

Please check with your healthcare professional for suitability.

Accessories - Wipes and Creams



Provox[©] Adhesive Remover

- Helps to remove adhesive comfortably and safely
- Reduces risk of skin irritation
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item N	o. Supplier	Description	Pack size
8012	Atos Medical AB	Provox Adhesive Remover	50



Provox[©] Cleaning Towel

- A convenient and hygienic way to clean around the stoma
- Small pack size, ideal for care bags, handbags and pockets
- No sting formula
- Does not leave any residue on the skin

ltem No.	Supplier	Description	Pack size
7244	Atos Medical AB	Provox Cleaning Towel	10



Provox[©] Skin Barrier

- Leaves a protective film on the skin to help prevent skin irritation
- Use prior to applying adhesive baseplate
- Convenient individual sachets
- No sting formula, safe for use around the stoma

Item No.	Supplier	Description	Pack size
8011	Atos Medical AB	Provox Skin Barrier	50



Spirotect[™] Skin Barrier

Item No.	Supplier	Description	Pack size
PTA500	Sheffmed Limited	Spirotect Skin Barrier	30



Spirotect[™] Adhesive Remover

Item No.	Supplier	Description	Pack size
PTA600	Sheffmed Limited	Spirotect Adhesive Remover	30



Cavilon© Skin Barrier

Item No.	Supplier	Description	Pack size
3343P	3M Health Care Ltd	Cavilon Skin Barrier Foam Applicator	5
3391G	3M Health Care Ltd	Cavilon Skin Barrier Cream 28g	1
3392G	3M Health Care Ltd	Cavilon Skin Barrier Cream 92g	1
3344E	3M Health Care Ltd	Cavilon Skin Barrier Wipe	30



Medihoney Barrier Cream

Item No.	Supplier	Description	Pack size
582	Derma Sciences	Medihoney barrier cream 50g	1

Accessories - Shower Aids



Provox[©] Life[™] Shower

For use with the Provox Life range.

Item No.	Supplier	Description	Pack size
8308	Atos Medical AB	Provox Life Shower	1



Provox[©] ShowerAid

For use with the Provox range.

Item No.	Supplier	Description	Pack size
7260	Atos Medical AB	Provox ShowerAid	1



Provox© Luna[™] ShowerAid

For use with Provox Luna Adhesive and Adhesive Strips.

ltem No.	Supplier	Description	Pack size
8016	Atos Medical AB	Provox Luna ShowerAid	1



Cascade Shower Protector

A flexible, waterproof cover which can be worn over the stoma.

ltem No.	Supplier	Description	Pack size
AS3835	Insight Medical	Cascade Shower Protector	2

Accessories - Foams Covers and Glue



Provox© **Tracheofix**

Item No.	Supplier	Description	Size	Pack size
1429H	Atos Medical AB	Provox Tracheofix 1 strip - Beige	55x60	10
1430H	Atos Medical AB	Provox Tracheofix 1 strip - Beige	70x70	10



ltem No.	Supplier	Description	Size	Pack size
1432H	Atos Medical AB	Provox Tracheofix 2 strips - Beige	55x60	10
1433H	Atos Medical AB	Provox Tracheofix 2 strips - Beige	70x70	10

Laryngofoam



Item No.	Supplier	Description	Size	Pack size
LALFF0023	Sheffmed Limited	Laryngofoam - Beige	63x65	30
LALFF0003	Sheffmed Limited	Laryngofoam - Beige	50x62	30

Purifoam

Item No.	Supplier	Description	Size	Pack size
CR3853	Insight Medical	Purifoam - Beige	65x65	30



Provox[©] Adhesive Strip

Protects Provox Life[™] Night Adhesive or Provox[®] Luna Adhesive when showering.

ltem No.	Supplier	Description	Pack size
8015	Atos Medical AB	Provox Adhesive Strip	10 pcs



Provox[©] Silicone Glue

ltem No.	Supplier	Description	Size	Pack size
7720	Atos Medical AB	Provox Silicone Glue	30ml	1

Accessories - Tape, Lubricating Jelly and Irrigation Fluid



Insil Tape

Item No.	Supplier	Description	Size	Pack size
INSIL200	Insight Medical	Insil Soft Silicone Perforated Tape	20mmx3m	1
INSIL400	Insight Medical	Insil Soft Silicone Perforated Tape	40mmx1.5m	1



Micropore Tape

ltem No.	Supplier	Description	Size	Pack size
153012	3M Health Care Ltd	Micropore tape	125mm x 5m	1



Optilube™

ltem No.	Supplier	Description	Size	Pack size
1120	Optimum Medical Solutions	Optilube Lubricating Jelly Sachets	5g	150 sachets
1121	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	42g	1
1122	Optimum Medical Solutions	Optilube Lubricating Jelly Tube	82g	1



Clinipod©

Clinipod is for cleaning purposes only and not the nebuliser.

ltem No.	Supplier	Description	Size	Pack size
CLINIPOD	Mayors Healthcare	Clinipod Sodium Chloride 0.9% Irrigation Solution	20ml	25



Irripod©

ltem No.	Supplier	Description	Size	Pack size
IRR24N	C D Medical Ltd	Irripod Sodium Chloride 0.9% Irrigation Solution	20ml	25

Normasol©



ltem No.	Supplier	Description	Size	Pack size
NOR206B	Molnlycke Health Care	Normasol Sterile Irrigation Fluid	25ml	20

Accessories - Dressing Packs



Sterile Dressing Pack Spec 10

ltem No.	Supplier	Description	Pack size
DP SPEC 10	Generic Supplier	Sterile Dressing Pack Spec 10	12 pcs



Sterile Dressing Pack Spec 35

Item No.	Supplier	Description	Pack size
DP SPEC 35	Generic Supplier	Sterile Dressing Pack Spec 35	12 pcs



Non-Woven Gauze Swabs

ltem No.	Supplier	Description	Pack size
SWAB SPEC 28	Generic Supplier	Non-Woven Gauze Swabs	100

Accessories - Gloves

Gloves



Item No.	Supplier	Description	Size	Pack size
GLO490B	Generic Supplier	Powder free nitrile gloves	Small	100
GLO489H	Generic Supplier	Powder free nitrile gloves	Medium	100
GLO488L	Generic Supplier	Powder free nitrile gloves	Large	100





Therabite® System Contains:

- 1 x Jaw Mobiliser
- 4 x Bite Pads
- 1 x Hand Aid
- 30 x Range Of Motion Scales
- 1 x Carry Bag
- 1 x Instructions For Use
- 1 x Patient Progress Log

TheraBite®

- For rehabilitation of reduced jaw opening (trismus).
- The TheraBite[®] Jaw Motion Rehabilitation System[™] is a clinically proven, anatomically correct and portable system. Designed to restore mobility after trismus and mandibular hypomobility.
- The TheraBite[®] Jaw Motion Rehabilitation System[™] can also help to reduce muscle pain by reducing joint inflammation. It can prevent reduced muscle strength during radiation treatment (chemotherapy) in the head and neck area.
- The TheraBite Jaw Motion Rehabilitation System is available in both adult and pediatric sizes.

ltem No.	Supplier	Description	Quantity
TH001	Atos Medical AB	TheraBite Jaw Motion Reh. System Adult	1
TH002	Atos Medical AB	TheraBite Jaw Motion Reh. System Pediat.	1
PA001	Atos Medical AB	TheraBite Bite Pad, Regular	4
PA002	Atos Medical AB	TheraBite Bite Pad, Edentulous	4
PA003	Atos Medical AB	TheraBite Bite Pad, Pediatric	4

Environmental policies



Fully recyclable packaging

To reduce our carbon footprint we use recyclable materials for every bit of our packaging.



Delivery using electric vehicles

Most of our deliveries are made using electric vehicles which is being expanded over the next 5 years.



Responsible forestry

Where possible, our paper and cardboard are sourced from FSC (Forest Stewardship Council) sources.



Biodegradable materials

Our packaging 'peanuts' are made from biodegradable corn starch and can be dissolved in water.

Contact us

There are several easy ways to get in touch:



Call or text us Tel: 0800 783 1659 Text: 0753 741 7928



Email us info@atos-care.co.uk



Visit our website www.atos-care.co.uk



Write to us Atos Medical, Cartwright House, Riverside Business Park, Tottle Road, Nottingham, NG2 1RT

Atos

Atos Care is the distribution arm of Atos Medical, the world-leading supplier of ostomy products for post laryngectomy and tracheostomy patients. Dedicated to making life easier for people living with a neck stoma, Atos Medical has decades of experience, and a deep understanding of what life is like after the operation.

So, although we supply the full range of products, we are much more than just a delivery service. Our expertise and focus mean that we can also provide full professional support and aftercare for you as an individual, wherever you are on your journey. We call this the Atos Circle of Care.

Atos Medical (UK), Cartwright House, Tottle Road, Nottingham, NG2 1RT Tel: 0800 783 1659 Text: 0753 741 7928 Email: info@atos-care.co.uk Web: www.atos-care.co.uk